


LE CUBE



 POLISHING


LE CUBE
Ref. 66250



LE CUBE

MANUAL POLISHING MACHINE

Simple, efficient and economical, Le CUBE is a manual polishing machine of Ø200 or 250mm, ideal for elementary polishing applications. With a reliable and robust construction, Le CUBE is the ideal partner for laboratories, schools and research centres.



CE

-  372 mm closed hood
700 mm opened hood
-  398 mm
-  420 mm
-  25 kg



Multijet watering ramp

- Adjustable flow
- Optimum lubrication of the support
- Evenly cools the platen

Motor with frequency inverter

- Powerful motor 180W
- High torque
- Silent functioning

Design & robustness

With sharp lines, this machine is easily recognizable with its futuristic, sleek and contemporary look.

Protective hood

- Transparent
- Easy to clean
- Protects from dust

Polishing platen

- Ø 200 or 250 mm
- 2 rotation directions (clockwise & anticlockwise)
- 2 rotation speeds (150 and 300 RPM)
- Accepts any type of polishing supports



Lighted aluminium buttons

buttons



TECHNICAL FEATURES

Removable polishing platen	✓
Motor power	180 W
Frequency inverter	✓
High-torque motor	✓
Speed and torque compensation	✓
Immediate stop of the platen (safety)	✓
Progressive acceleration	✓
Direct transmission on bearings	✓
Anti-splash flange	✓
Voltage / Frequency	230 V/50-60 Hz